


*Crusade*TM



KANEKA
MEDICAL
PRODUCTS

Crusade™

New possibility in the era of DES for bifurcation lesion.


1 Optimized configuration and material achieved high shaft maneuverability.

- Thin and flexible tip : **Distal shaft**
- Pliable and sturdy : **Proximal shaft**

2 Configuration of double-layered lumens allows more flexible GW movement

Product specifications

Catalogue Number	Maximum Guide Wire Outer Diameter	Catheter Effective Length	
KMF0114A	0.014inch (0.36mm)	1400mm	



[Precautions during use]

Inserting a highly penetrating guidewire while this catheter is left in a tortuous lesion may break the catheter shaft resulting in blood vessel injury.

⊗ Do not re-use

STERILE EO

Sterilized using ethylene oxide



Consult instructions for use

Manufacturer

KANEKA CORPORATION

3-18, 2-Chome, Nakanoshima, Kita-ku, Osaka-city, OSAKA, 530-8288 JAPAN

Tel.No. : (+81)-(0)6-6226-5256 **CONTACT** Tel.No. : (+81)-(0)3-5574-8136

EC-Representative

KANEKA PHARMA EUROPE N.V.

Nijverheidsstraat 16, 2260 Westerlo-Oevel, Belgium

Tel.No. : (+32)-(0)14-256-297 Fax No. : (+32)-(0)14-256-298

CE 0123

Ref. No.E107-59e Date Mar. 2016 2000 IT