




# Thrombuster™ II

6Fr  
7Fr

# Thrombuster™ II

## Features



## Product specifications

Catalogue Number	Guide Catheter Minimum ID Compatibility	Maximum Guide Wire Outer Diameter	Shaft Diameter				Guidewire Lumen Part		Catheter Effective Length
			Outer		Inner		Maximum Outer Diameter	Length	
			Distal	Proximal	Distal	Proximal			
T2R6SA4W	0.070inch (1.78mm)	0.014inch (0.36mm)	3.9Fr (1.30mm)	3.9Fr (1.30mm)	3.0Fr (1.00mm)	3.3Fr (1.10mm)	5.1Fr (1.70mm)	10mm	1400mm
T2R7SA4W	0.080inch (2.03mm)	0.014inch (0.36mm)	4.5Fr (1.50mm)	4.6Fr (1.52mm)	3.6Fr (1.20mm)	3.9Fr (1.32mm)	5.7Fr (1.90mm)	10mm	1400mm





Do not re-use



Sterilized using ethylene oxide



Consult instructions for use



Manufacturer

**KANEKA CORPORATION**

3-18, 2-Chome, Nakanoshima, Kita-ku, Osaka-city, OSAKA, 530-8288 JAPAN

Tel.No. : (+81)-(0)6-6226-5256 [CONTACT](#) Tel.No. : (+81)-(0)3-5574-8136

EC-Representative

**KANEKA PHARMA EUROPE N.V.**

Nijverheidsstraat 16, 2260 Westerlo-Oevel, Belgium

Tel.No. : (+32)-(0)14-256-297 Fax No. : (+32)-(0)14-256-298

